




**+ Pekkton® ivory.**

high performance polymer for  
definitive aesthetic restorations  
on implants.

**anaxDENT**


# High performance polymer PEKK.

The material of the future.



## The performance pyramid

PEKK stands at the apex of the PAEK family.



## The PEKK chemical structure

The PEKK chemical structure guarantees the best mechanical properties of all PAEK. The PEKK has for example up to 80% greater compression strength than PEEK.

## The latest generation: PEKK

The superior high performance polymer Pekkton® ivory is the solution for definitive, aesthetic, patient-friendly restorations in implant dentistry.

To meet the demanding biomechanical criteria of prosthetic solutions, Cendres+Métaux has secured the exclusive rights to Pekkton® (PEKK): the top product among thermoplastic polyaryletherketones.

High performance polymer:



- Has similar characteristics to those of human bone
- Shock absorption
- Light-weight and natural for the patient

Pekkton® ivory:

- Available exclusively from Cendres+Métaux
- The latest generation of polymers

# Pekkton® ivory.

Similar properties to biological materials.



## Biological

Pekkton® ivory has similar characteristics to those of human bone.

## Compressive strength

Pekkton® ivory is comparable to human dentine and cortical bone in terms of compressive strength. It allows better bio-mechanical integration than typical non-precious metals.

## Unique properties

Thanks to its unique, nature-adapted, mechanical, physical and chemical properties, Pekkton® ivory offers a broader spectrum of applications than other polymers.

Compressive strength	246	MPa
Bending strength	200	MPa
Flexural modulus	5.0	GPa
E-modulus	5.1	GPa
Tensile strength (MPa)	115	MPa

Melting point	363	°C
Density	1.4	g/cm <sup>3</sup>
Water absorption	8.7	µg/mm <sup>3</sup>
Solubility	0.2	µg/mm <sup>3</sup>
Hardness (DIN EN ISO 2039-1)	252	MPa

### Country specific indication may apply.

The current FDA indication applies specifically to the pressable version of Pekkton® as a framework material for long-term, temporary crowns and bridges. The manufacturer is currently pursuing further FDA clearance to align with current CE (European) indications for use as a final prosthesis in both pressable AND millable formats.

# Pekkton® ivory.

For successful restorations on implants.



## Modern materials and traditional techniques

Important aspects for this case on 6 implants were the feel of the restorations, the aesthetics, youthful appearance, depth of colour and overall symmetry. Pekkton® ivory was the patients' material of choice.

- a) Pressed Pekkton® ivory framework, prepared for veneering
- b) Individual pressed copings veneered with ceramic
- c) Natural gingival design with light-cured composite





### Full-arch bridge on implants

Pekkton® ivory framework provides the patient with a lightweight aesthetic and extremely strong restoration that offers shock absorbing properties not found in traditional implant restorative materials.

### Other indications



Removable restorations, such as secondary constructions on bars and telescopic crowns, transversal connectors, occlusal splints and denture bases.



Definitive supported, veneered single crowns and bridges with a maximum of one pontic on natural teeth.




Unveneered parts e.g. crown margins and backings.

## Clinical cases.

Results in day-to-day life.

### Full-arch bridge screwed on 5 implants



### Molar crown on an implant



**Skeletal partial framework stabilised by 5 natural telescopic abutments in the upper jaw**



**Screwed restoration on implants in the upper and lower jaw**



# Pekkton® ivory.

## Benefits.

### Shock-absorbing

Excellent for your restoration on implants



### For milling and pressing

Easy to mill and perfect to polish



### Break new ground, join the new world of materials with Pekkton® ivory

- Shock-absorbing
- Metal-free
- Similar properties to biological materials
- Extremely low water absorption
- Caries control
- Visible in the X-ray
- Excellent «weight-to-strength ratio»
- Can be sterilised
- Biocompatible

### Pekkton® ivory: analog or digital processing

- For an excellent fit and full control in the dental laboratory
- Time savings and guaranteed consistent quality throughout the digital process



## Aesthetic

High flexibility in the choice of veneering



## Definitive

Happy patients thanks to the natural feel in the mouth



### A wide range of veneering

- Composites
- Shells
- Teeth
- Ceramic crowns

### Pekkton® ivory makes your patients smile

- Strong and light
- No metal allergies
- No taste
- No thermal or electrical conductivity

# Traditional processing. Pressing.

For an excellent fit and full control in the dental laboratory



## PEKKtherm

The device makes it easy and safe to stabilise the muffle temperature (out of the preheating furnace) to the correct pressing temperature. Pekkton® ivory is then melted before the pressing process.

## PEKKpress

This device is used to efficiently and gently inject Pekkton® ivory after the material and muffle have been prepared in the PEKKtherm device.



PEKKtherm



PEKKpress

# Digital processing.









## Milling.

Time savings and guaranteed consistent quality throughout the digital process



# Portfolio.

The essence of success.

Order No. Cendres + Métaux		Description
01060011		Pekkton® ivory milling blank ø 98.5/16mm (with steps)
01060020		Pekkton® ivory milling blank ø 98.5/20mm (with steps)
01060022		Pekkton® ivory milling blank ø 98.5/24mm (with steps)
01060028		Pekkton® ivory milling blank ø 95/16mm (compatible with Zirkozahn®)
01060030		Pekkton® ivory milling blank ø 95/20mm (compatible with Zirkozahn®)
01060032		Pekkton® ivory milling blank ø 95/24mm (compatible with Zirkozahn®)
01060003		Pekkton® ivory – press ingots / 10 p.
70202332		Disposable press-stamp (ø 12mm, incl. seal) / 50 p.
08000626		Disposable press-stamp (ø 12mm) / 50 p.
08000627		Disposable press-stamp (ø 26mm) / 20 p.
08000628		PEKKpress muffle former set 200g
08000629		PEKKpress muffle former set 600g
70200086		AnaxVest PM 3kg (30 x 100g) incl. Liquid 1L
70100076		Anaxdent – Bond LC 10ml
70202393		PEKKpress – pressing unit
70202394		PEKKtherm – temperature stabilisation and melting furnace

Country specific indication may apply.  
Please follow the instructions for use for further  
information on indications, contraindications,  
warnings and precautions.

Pekkton® ivory is based on OXPEKK® by OPM,  
Oxford Performance Materials, Inc., USA.

## anaxDENT

**CENDRES<sup>+</sup>  
MÉTAUX**

anaxdent North America  
1408 W. Broadway St.  
Ardmore, OK 73401, USA

Cendres + Métaux SA  
Rue de Boujean 122  
CH-2501 Biel/Bienne



[www.facebook.com/dentcmsa](http://www.facebook.com/dentcmsa)



[www.youtube.com/user/dentalcmsa](http://www.youtube.com/user/dentalcmsa)



[www.pekkton.com](http://www.pekkton.com)

Phone +1 877 897 65 98  
[usinfo@anaxdent.com](mailto:usinfo@anaxdent.com)

[www.anaxdentusa.com](http://www.anaxdentusa.com)

Phone +41 58 360 20 00  
Fax +41 58 360 20 15  
[info@cmsa.ch](mailto:info@cmsa.ch)

[www.cmsa.ch/dental](http://www.cmsa.ch/dental)